

Optional: Spreizersystem für die lumbale Bandscheibe
 Optional: Spreading system for the lumbar intervertebral disc
 Opcional: Sistema de separadores para el disco intervertebral lumbar
 Optionnel: Système écarteur pour le disque intervertébrale lombaire
 Opzionale: Sistema di divaricazione per il disco intervertebrale lombare



23.20.02



57.58.35
57.58.36
57.58.37
57.58.38
57.58.39
57.58.40
57.58.42



57.57.52



57.60.63
57.60.64
57.60.65
57.60.66
57.60.67
57.60.68
57.60.69



57.60.83
57.60.84
57.60.85
57.60.86
57.60.87
57.60.88
57.60.89



68.18.93



49.50.10



49.97.02



49.97.90

PART 5

Item number	Description	Quantity
23.20.02	KOCHER Retractor · 12 × 45 mm · 21 cm	1
57.58.35	Speculum for counter-retractor · 40 mm	1
57.58.36	Speculum for counter-retractor · 45 mm	1
57.58.37	Speculum for counter-retractor · 50 mm	1
57.58.38	Speculum for counter-retractor · 55 mm	1
57.58.39	Speculum for counter-retractor · 60 mm	1
57.58.40	Speculum for counter-retractor · 65 mm	1
57.58.42	Speculum for counter-retractor · 75 mm	1
57.57.52	Counter retractor with joint	1
→ 57.60.63	PEEK Blade · semi transparent to X-ray · 12 mm · 40 mm	1
→ 57.60.64	PEEK Blade · semi transparent to X-ray · 12 mm · 45 mm	1
→ 57.60.65	PEEK Blade · semi transparent to X-ray · 12 mm · 50 mm	1
→ 57.60.66	PEEK Blade · semi transparent to X-ray · 12 mm · 55 mm	1
→ 57.60.67	PEEK Blade · semi transparent to X-ray · 12 mm · 60 mm	1
→ 57.60.68	PEEK Blade · semi transparent to X-ray · 12 mm · 65 mm	1
→ 57.60.69	PEEK Blade · semi transparent to X-ray · 12 mm · 75 mm	1
57.60.83	PEEK Blade · excellent X-ray transparency · 8/12 mm · 40 mm	1
57.60.84	PEEK Blade · excellent X-ray transparency · 8/12 mm · 45 mm	1
57.60.85	PEEK Blade · excellent X-ray transparency · 8/12 mm · 50 mm	1
57.60.86	PEEK Blade · excellent X-ray transparency · 8/12 mm · 55 mm	1
57.60.87	PEEK Blade · excellent X-ray transparency · 8/12 mm · 60 mm	1
57.60.88	PEEK Blade · excellent X-ray transparency · 8/12 mm · 65 mm	1
57.60.89	PEEK Blade · excellent X-ray transparency · 8/12 mm · 75 mm	1
68.18.93	DOCKHORN Retractor · 55 × 14 mm · 21 cm	1
49.50.10	Mediccontainer · silver · 290 × 290 × 135 mm	1
49.97.02	Wire basket · 240 × 255 × 50 mm	1
49.97.90	Wire basket · 240 × 250 × 40 mm	1